**PROCTOSCOPE
WITH LIGHT GUIDE**

Product Description and Application

Description	Model Size	Product Code
PROCTOLUX®	19mm	702000
PROCTOLUX®	25mm	702010

PROCTOLUX™ is a range of single use, fully disposable and self-illuminating rigid endoscope. These devices are used during either rectal examination or rectal surgery for use in examination, diagnosis and treatment in anal cavity, rectum, or sigmoid colon. The PROCTOLUX™ Proctoscopes come in two sizes 19mm and 25mm diameters. All PROCTOLUX™ devices are shipped ready for use.

Should light output become impaired in use due to debris or smoke obscuring the light output the end of the light guide can be cleared by wiping.

Re-use of this device represents an unnecessary and avoidable risk of cross infection to the patient.

Application

The G&N PROCTOLUX™ Proctoscope and associated Illuminator are intended to provide a means of rigid speculum to the anal cavity, rectum, or sigmoid colon and to supply effective illumination thus permitting direct visual examination and treatment by a suitably qualified and experienced healthcare professional.

The devices key functions are:

1. To enable trained health and care professionals, to effectively view anal cavity, rectum, or sigmoid colon to perform invasive investigational procedures.
2. Reduce the complexity and increase the ease of use by having a light source directly attached to the device enabling the area to be viewed by the operator facilitated by an internal light source that can illuminate the area without the need for extra equipment.

Intended Use

The G+N Medical 19mm and 25mm PROCTOLUX™ Proctoscopes are intended to provide means of rigid endoscopic access to the patient's rectum and supply effective illumination to permit direct visual examination by a suitably qualified and experienced healthcare professional.

Indication and Patient Group

The Griffiths and Nielsen 19mm and 25mm PROCTOLUX™ Proctoscopes are indicated for use when a suitably qualified and experienced healthcare professional requires to perform either a direct rectal examination or surgery of the non-contraindicated patient using a suitable rigid endoscope.

Instructions for Use

1. When ready to use the device, remove the device from its packaging Separate the peel pouch and remove the product. Discard the package. Examine the pack before use. The pack shall not be issued nor the proctoscope used if the wrapper or seal is broken DO NOT USE IF PACKAGING OR PROCTOSCOPE IS DAMAGED
2. To activate the light, pull out the blue tab in the light pack in the device handle. Pull out the blue tab in the light pack in the device handle to activate the light source before insertion into the patient. Maximum illumination intensity is maintained for a minimum of 30 minutes. The light pack will operate for a considerably longer period, with gradually diminishing intensity. Once activated the light cannot be switched off.
3. Ensure that the appropriate size of proctoscope is selected before commencing the procedure
4. Apply a suitable lubricant to the outer aspect of the proctoscope and then gradually and gently insert the proctoscope into the anal canal and rectum
5. When the proctoscope is in the required position, remove the inner obturator to facilitate insufflation (if required) and observation of the rectum
6. Upon completion of the examination, gradually and gently withdraw the proctoscope from the anal canal and then confirm that the patient is comfortable
7. After use, ensure that the product is disposed of safely into an appropriate clinical waste container in accordance with local policies and procedures as described below.
8. Dispose the device in an appropriate way as described in the disposal section below.

IMPORTANT Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Contraindications

Acute peritonitis	Acute Severe Diverticulitis	Diverticular Abscess
Colonic Necrosis	Recent Colonic Surgery	Anal Fissure
Fulminant Colitis	Toxic Megacolon	

	Do not re-use		Do not use if package is damaged
	Catalogue number		Keep dry
	Batch number		Shock protection (Type BF Applied Part)
	Use by date		Federal or EU Law restricts this device to be sold by or on the order of a licensed Health Practitioner
	Date of manufacture		Warning/caution consult accompanying documents
	Direct electrical current at voltage shown		Consult instructions for use/consult operating instructions
	Not made with natural rubber latex		Non-sterile
	Do not look directly into the light source		Manufacturer and address of the manufacturer
	Fragile – handle package carefully		Number of items per pack
	No DEHP		IP rating
	Medical Device		Medical Device in conformity with Regulation (EU) 2017/745

GMDN 35255 Global Medical Device Nomenclature

Disposal

Discard PROCTOLUX™ with the light pack still attached, as “clinical waste” according to hospital, or local protocols. Note: If required the light pack can be removed by firmly pulling the entire pack downwards. Do be careful as it may well release suddenly. Dispose of the light pack responsibly.

Warnings

Before each use carefully examine the instrument to ensure that there are no visible scratches, chips, voids, etc. These could cause the instrument to lose its stability potentially compromising safety. The patient should be informed when the Proctoscope is to be introduced and removed.

Adjustments may be necessary to provide improved comfort.

WARNING		Inspect each device and packaging for damage prior to use
WARNING		For hygienic reasons always wear surgical gloves when handling the device before and after use. The device will come into contact with biological fluids which are potentially infectious.
WARNING		When in transit or storage, device may be subject to damage beyond the control of the manufacturer or supplier.
WARNING		Never use device with laser equipment.
WARNING		Treat used devices as biohazardous infectious material. Dispose of used devices in suitable disposal unit or in accordance with local regulations.
WARNING		Device is not compatible with any sterilization method.
WARNING		No modification of this equipment is allowed.
WARNING		Never use an ultrasonic cleaner for the device, the device should not be cleaned or sterilised in any way as it provided in clean hygienic form for single use only.

FOR SINGLE USE ONLY:

The G&N Medical device is designed for single use only; do not reuse device. G&N does not have data regarding reuse of this device. Reuse may cause device failure or procedural complications including device damage, compromised device biocompatibility, and device contamination. Reuse may result in infection, serious injury, or patient death. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.

Cautions:

For use by trained healthcare personnel only

Federal (USA), UK and European Union Law restricts this device to sale, distribution and use by or on the order of a physician.

Storage

Store devices in their original packaging at room temperature and normal air humidity. The shelf life is indicated on the product.

CS Lifesciences Europe Limited
3 Inns Quay, Dublin 7, Ireland

110002/3
05/2021

Griffiths and Nielsen Ltd.
Maydwell Avenue, Stane Street,
Slinfold, Horsham, RH13 0GN, UK

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