




FITLEGS[®]2



ASK FOR COMFORT.

ASK FOR FITLEGS®2

G+N FITLEGS®2 compression stockings are designed to provide pressure to the legs to help improve circulation and drainage. Featuring our DUOFIT™ technology, the unique construction of the weave allows us to vastly reduce the amount of elastane, whilst still providing the amount of compression required. This means that FITLEGS®2 are softer to use, making them more comfortable to wear over long periods of time.



DUOFIT™ technology provides a two-way stretch for easier application.



Available in below knee and thigh length and tights.



Ankle and calf measurement for Below Knee along with Thigh circumference measurement for tights and Thigh Length ensuring correct fit for patient.



Colour coded packaging reduces time identifying correct size.



Four sizes within each length reduces inventory.

WHY USE FITLEGS 2?

FITLEGS 2 Compression garments have a compression value of 23-32mmHg.

- An essential part of ongoing care and should be worn regularly to maintain reduction in limb volume
- Specific compression level allows for venous return and lymphatic drainage
- FITLEGS 2 Below knee and Thigh length available in Black, White and Beige. Tights available in Black and Beige.



INTELLIGENT SOLUTIONS FOR LYMPHODEMA MANAGEMENT

With only 4 sizes within each style, G&N FITLEGS 2 allow hospitals and clinics to benefit from a reduced inventory – saving space in stockrooms and Healthcare Professionals time when identifying the correct size.

Only 4 sizes
with colour coded
packaging

PRODUCT CODES & SIZING GUIDE





FITLEGS® 2 (23-32mmHg)

Suitability Check

	Ankle measurement	Size	BELOW KNEE Order Code	THIGH Order Code	Measurement based on height & weight	TIGHTS Order Code	Maximum Thigh Circumference
Knee	19 - 23cm	Small	457030	457034		457046	65cm
	23 - 26cm	Medium	457031	457035	See table	457047	70cm
	26 - 29cm	Large	457032	457036	below for	457048	78cm
	29 - 33cm	Ex-Large	457033	457037	full sizing for	457049	82cm
Below knee	19 - 23cm	Small	457016	457021	Fitlegs Class	457050	65cm
	23 - 26cm	Medium	457017	457022	2 tights	457051	70cm
	26 - 29cm	Large	457018	457023	range	457052	78cm
	29 - 33cm	Ex-Large	457019	457024		457053	82cm
White	19 - 23cm	Small	457038	457042			65cm
	23 - 26cm	Medium	457039	457043			70cm
	26 - 29cm	Large	457039	457044			78cm
	29 - 33cm	Ex-Large	457040	457045			82cm

Wearing the correct size compression stockings is important and customers must be measured by a clinician before purchasing our stockings.



Griffiths and Nielsen Ltd
 Maydwell Avenue, off Stane Street,
 Slinfold, Horsham, West Sussex RH13 0GN
 Tel: +44 (0) 1403 799 190
 Fax: +44 (0) 845 263 8907
 medsales@gandn.com

www.gandn.com

