

Your patients deserve it all. Get the **NEXADVANTAGE**.

**Experience the Most Positive Thing to Happen to Negative Pressure Wound Therapy** 





# What is the **NEXA**DVANTAGE?

### The NEXA Negative Pressure Wound Therapy (NPWT) System

has been developed to provide uncompromising negative pressure therapy for patients with chronic and acute wounds.

### SIMPLE:

- With no complicated interfaces or touchscreens, NEXA's ultra-simple, one-button operation requires minimal training and reduces the potential for user error
- Sleek, streamlined design with no maintenance required
- Accurate pressure control is built in, as the pre-set SmartFlow<sup>™</sup> tubing connector constantly regulates negative pressure

### **QUIET:**

- Pump runs smoothly and consistently
- No audible notifications means less "alarm anxiety" for staff, patients and family

### PORTABLE:

- Lightweight and small, ideal for mobile patients
- Easy to change, compact, lightweight pouch that collects fluid away from the wound
- Handy carrying case for patients on the go
- Battery lasts 10 hours before recharging

### AFFORDABLE:

- Designed to operate for up to 8 weeks (56 days or 1344 hours)
- Usable hours are depleted only when the device is providing therapy
- Remaining hours can be used at your convenience on the next patient
- Immediate availability with this off-the-shelf device that requires no tracking and incurs no runaway rental costs
- Disposable device eliminates servicing and maintenance costs



### **Smart Choices Add Up**

When you're focused on effective, efficient patient care, the right partner can make a huge difference. NEXA is the affordable NPWT option for your facility.

#### Here's why:



\*For device only. Additional costs for dressings and EcoPouch™ containers will vary depending on the number required for each patient.

A Powerful Balance of SIMPLICITY and PERFORMANCE

This multi-week, multi-patient, disposable system has three innovations that make all the difference.



### Hydrodynamic Fluid Pump

Continues to remove fluid even in the presence of a minor dressing air leak

- Constantly, actively removes exudate from wound, reducing potential for skin maceration
- No risk of internal device contamination because the disposable pump never touches the inside of the tubing or fluid container pouch



#### EcoPouch™

Draws exudate away from the wound site into the separate, disposable, flexible pouch

- 120 ml capacity
- Much less shelf space required
- No sterilisation needed
- Requires 58% less polymer content than conventional rigid canisters<sup>†</sup>



### **SmartFlow™ Connector**

Pressure is constantly maintained at the wound

- Maintains negative pressure levels preset to -125mmHg
- No need for complex pressure monitoring sensors or feedback loops within device, dressing, or tubing



## **ORDERING** Information

Product	Description	Part#	UOM
NEXA NPWT System Device Kit	NEXA therapy unit, power supply, carry case	400-00-300-GBR	1
NEXA Dressing -125mmHg	1 occlusive drape, 1 tubing connection pad, 1 foam filler (15cm x 7.5cm)	414-00-501	Box of 5
NEXA EcoPouch™ Fluid Container Pack	120ml flexible fluid container pouch and tubing	414-00-502	Box of 5
NEXA Drape	Occlusive drape individual sheets	414-00-600	Box of 10
Carrying Case	Carrying case for NEXA therapy unit	414-00-400	1
Power Supply	Power supply for NEXA therapy unit	414-00-451	1





#### S&N

Maydwell Avenue, Stane Street, Slinfold, Horsham, West Sussex RH13 OGN

Tel: +44 (0) 1403 799190 Email: sales@gandn.com

www.gandn.com



