


Advanced Wound Care & Burn Dressings

The Silverlon® Advantage

The Silverlon Antimicrobial Silver-Plated Dressing technology provides the highest amount of silver ions for immediate antimicrobial activity to reduce the risk of infection.



Silverlon® dressings provide visible results.

Performance01

Power02

Speed03

With Silverlon®, you can reference multiple robust, peer-reviewed clinical studies across multiple specialties.

Silverlon® Wound Dressings contain more metallic silver than any other dressing.

Silverlon® produces a broad-spectrum antimicrobial effect that rapidly kills microorganisms, eradicating MRSA in as little as 4 hours and keeps working for up to 7 days.

What makes Silverlon® dressings unique?

- Sustained release01
- Deliverable silver ions02
- No staining03

Silverlon® silver plated nylon technology offers many unique advantages over other silver dressings.




Silverlon® SurfaceSilverlon® Magnified

Silverlon® dressings contain more metallic silver than other silver dressings

Proprietary Name	Ag Content mg/100 cm²
SILVERLON®	546
Askina® Calgitrol® Ag+	141
ACTICOAT◊	104
Contreet Foam	85
Contreet	32
AQUACEL® Ag	8.3
SilvaSorb®	5.3
ACTISORB™ Silver 220	2.7

ACTICOAT◊ is a registered trademark of Smith+Nephew
Askina® Calgitrol® Ag+ is a registered trademark of B. Braun
AQUACEL® Ag is a registered trademark of ConvaTec Inc.
SilvaSorb® is a registered trademark of Virchow Biotech, Inc.
ACTISORB™ is a trademark of Medical Action Industries, Inc.

Silverlon® has an antimicrobial effect within 2 hours against MRSA in the dressing, with eradication in 4 hours.



Time and CFU/ml		0.0Hr	0.5Hr	1.0Hr	2.0Hr	4.0Hr	6.0Hr	24.0Hr
Bacterial species								
Staphylococcus aureus	MRSA ATCC 33591	8,250,000	2,600,000	270,000	500	0	0	0
Pseudomonas aeruginosa	ATCC 9027	9,000,000	4,000,000	0	0	0	0	0

Other “Silver” dressings

Silver impregnated dressings

- May have as little as 1/50th as much silver as Silverlon®
- Silver layer may not be in direct contact with the wounded tissue.
- Typically supported by in vitro data only.

Silver nanocrystalline dressings

- Typically 1/5th the silver of Silverlon®
- Staining and stinging can occur

Indications

Silverlon® Wound Contact Dressings, Burn Contact Dressings, and Elastic Burn Wraps are professional use wound dressings indicated for local management of:

- 1st and 2nd degree burns
- Donor and graft site wounds
- Traumatic and surgical wounds
- Diabetic, pressure and venous ulcers
- Infected wounds*

*Silverlon® Wound Contact Dressings, Burn Contact Dressings, and Elastic Burn Wraps may be used on those infected wound types listed above. When used on infected wounds, the infection should be treated as per local clinical protocol.

Silverlon® Wound Contact Dressings, Burn Contact Dressings, and Elastic Burn Wraps are indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.

Contraindications

- Avoid using Silverlon® Wound Contact Dressings, Burn Contact Dressings, and Elastic Burn Wraps on patients with known sensitivity to silver or nylon.
- Silverlon® Wound Contact Dressings, Burn Contact Dressings, and Elastic Burn Wraps are not intended for use on 3rd degree burns.

Our Silverlon® products

Dressings are available in multiple shapes & sizes

Silverlon® Flex Dressings



Product and Catalogue No.	Size
WCD44-NPS-FLX	10 cm x 10 cm
WCD48-NPS-FLX	10 cm x 20 cm
WCD88-NPS-FLX	20 cm x 20 cm
WCD22-NPS-FLX	5 cm x 5 cm
WCD22-NPS	6 cm x 6 cm
WCD45-NPS	10 cm x 13 cm
WCD58-NPS	13 cm x 20 cm
WCD1212-NPS	30 cm x 30 cm
WCD88-NPS	20 cm x 20 cm

Silverlon® Negative Pressure Dressings



Product and Catalogue No.	Size
NPD-22	5 cm x 5 cm
NPD-45	10 cm x 13 cm
NPD-58	13 cm x 20 cm
NPD-512	13 cm x 30 cm
NPD-1212	30 cm x 30 cm

Silverlon® Burn Wraps



Product and Catalogue No.	Size
BWD-466	10 cm x 168 cm
BWD-6108	15 cm x 274 cm

Silverlon® Burn Chest Dressings




Product and Catalogue No.	Size
BCD-CDL	76 cm x 65 cm
BCD-CDM	66 cm x 64 cm

Silverlon® Burn Contact Dressings



Product and Catalogue No.	Size
BCD-44	10 cm x 12 cm
BCD-48	10 cm x 20 cm
BCD-816	20 cm x 40 cm
BCD-1616	41 cm x 41 cm
BCD-2424	61 cm x 61 cm

Silverlon® Burn Pad Dressings



Product and Catalogue No.	Size
BPD-44	10 cm x 12 cm
BPD-48	10 cm x 20 cm
BPD-816	20 cm x 40 cm

Silverlon® Wound Pad Dressings



Product and Catalogue No.	Size
WPD-22	5 cm x 5 cm
WPD-23	5 cm x 8 cm
WPD-26	5 cm x 15 cm
WPD-28	5 cm x 20 cm
WPD-210	5 cm x 25 cm
WPD-212	5 cm x 30 cm
WPD-38	8 cm x 20 cm
WPD-310	8 cm x 25 cm
WPD-316	8 cm x 40 cm
WPD-44	10 cm x 12 cm

Silverlon® Wound Packing Strips



Product and Catalogue No.	Size
WPS-112	3 cm x 30 cm
WPS-124	3 cm x 61 cm

Silverlon® Surgical Island Dressings



Product and Catalogue No.	Size
ID-23	5 cm x 8 cm
ID-34	8 cm x 10 cm
ID-44	10 cm x 10 cm
ID-46	10 cm x 15 cm
ID-410	10 cm x 25 cm
ID-412	10 cm x 30 cm
ID-414	10 cm x 36 cm
ID-66	15 cm x 15 cm

Silverlon® Wound Contact Dressings



Product and Catalogue No.	Size
WCD-22	5 cm x 5 cm
WCD-44	10 cm x 12 cm
WCD-412	10 cm x 30 cm
WCD-1012	25 cm x 30 cm
WCD-466	10 cm x 168 cm

Silverlon® Catheter Dressings



2.5 cm Disc Dressings	Hole Size
IVCD-15	1.5 mm
IVCD-40	4 mm
IVCD-70	7 mm



Sold & distributed by:

G&N

Maydwell Avenue, Stane Street,
Slinfold, Horsham, West Sussex RH13 OGN

Tel: +44 (0) 1403 799190

Email: sales@gandn.com

www.gandn.com

argentum
medical

A Bravida Medical Company