

Lifesaver reference guide

Silverlon® Lifesaver Catheter Dressings

IVCD-15

- 2.5 cm Dressing
- 1.5 mm I.D.
- Up to 4.5 French
- Ideal for Central Lines



IVCD-40

- 2.5 cm Dressing
- 4.0 mm I.D.
- Up to 12 French
- Ideal for Central Lines, Ortho Fixator Pins



IVCD-70

- 2.5 cm Dressing
- 7.0 mm I.D.
- Up to 21 French
- Ideal for Dialysis, Large Multi Lumen Catheters and Large MAC (Multi Access Catheters)



Applications

Vascular Applications

- CVC/PICC Lines
- Huber Needles
- Dialysis Catheters
- Arterial Catheters
- Epidural Catheters
- Peripheral I.V. Catheters
- Mechanical Circulatory Support
- Broviacs

Specialty Applications

- External Fixator Pins (LLRS, Podiatry, etc)
- Chest Tubes
- Gastrostomy Feeding Tubes
- Non-Vascular Percutaneous Devices
- Burns
- Chest Tubes
- Non-Intact Skin
- Chest Tubes

Indications

Silverlon® Lifesaver IV/ Catheter Dressings are professional use wound dressings indicated for local management of:

- Vascular access sites (including around IV catheters, central venous lines, arterial catheters, dialysis catheters);
- Drainage tube access sites;
- Epidural catheter access sites;
- Infected wounds*

Contraindications

Avoid using Silverlon® Lifesaver IV/Catheter Dressings on patients with known sensitivity to silver or nylon.
For complete Instructions for Use, please visit www.silverlon.com.

* Silverlon® Lifesaver IV/Catheter Dressings may be used on those infected wound types listed above. Where used on infected wounds, the infection should be treated as per local clinical protocol.
Silverlon® Lifesaver IV/ Catheter Dressings are indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.



Ordering Information

Product and Catalogue No.	Size	Qty
Silverlon® Catheter Dressings		
IVCD-15	IV/Cath Dressings - 2.5 cm Lifesavers 2.5 cm Dia 1.5 mm access	10
IVCD-40	IV/Cath Dressings - 2.5 cm Lifesavers 2.5 cm Dia 4.0 mm access	10
IVCD-70	IV/Cath Dressings - 2.5 cm Lifesavers 2.5 cm Dia 7.0 mm access	10

To learn more about Silverlon® solutions, please contact:



G&N
Maydwell Avenue, Stane Street,
Slinfold, Horsham, West Sussex RH13 OGN
Tel: +44 (0) 1403 799190
Email: sales@gandn.com
www.gandn.com