






# SUPERIOR EASE OF USE AND COMFORT

G+N White anti-embolism stockings help to prevent the formation of Deep Vein Thrombosis (DVT) in the legs by applying clinically proven graduated compression.<sup>1,2</sup> G+N White anti-embolism stockings are designed to be simple to fit and apply, which reduces sizing errors and nursing time.



DUOFIT™ technology provides a two-way stretch for easier application.



Improved inspection hole design and placement for superior comfort.



Single ankle measurement reduces nursing time when measuring patients.



Available in below knee and thigh length.



Six sizes within each length reduces inventory.





Anti-bacterial fabric for odour management.




Colour coded packaging reduces time identifying correct size.

## HOW DOES IT WORK?

The compression within the garment is highest at the ankle and reduces further up the leg. Pressure forces blood back towards the heart.



With G+N White anti-embolism stockings



Without G+N White anti-embolism stockings



## INTELLIGENT SOLUTIONS FOR DVT PREVENTION

With only 6 sizes within each style, G+N White anti-embolism stockings allow hospitals to benefit from a reduced inventory – saving space in stockrooms and nursing time when identifying the correct size.

**Only 6 sizes**  
with colour coded  
packaging

1. Independent Sigel Certification. Data on File.  
2. Sigel B., et al. Type of Compression for Reducing Venous Stasis. Archives of Surgery. February 1975. Vol 110: 171-175.

# PRODUCT CODES & SIZING GUIDE



## FITLEGS® AES White


Suitability Check

Ankle measurement	Size	BELOW KNEE Order Code		THIGH Order Code		Maximum Thigh Circumference
16 - 19cm	Slimline		403004			
19 - 23cm	Small		404000		404004	65cm
23 - 26cm	Medium		404001		404005	70cm
26 - 29cm	Large		404002		404006	78cm
29 - 33cm	Ex-Large		404003		404007	82cm
33 - 35cm	Bariatric		403006			

## FITLEGS® AES GRIP White

Suitability Check

Ankle measurement	Size	BELOW KNEE Order Code		THIGH Order Code		Maximum Thigh Circumference
16 - 19cm	Slimline		403005			
19 - 23cm	Small		404008		404012	65cm
23 - 26cm	Medium		404009		404013	70cm
26 - 29cm	Large		404010		404014	78cm
29 - 33cm	Ex-Large		404011		404015	82cm
33 - 35cm	Bariatric		403007			



### Griffiths and Nielsen Ltd

Maydwell Avenue, off Stane Street,  
Slinfold, Horsham, West Sussex RH13 0GN

Tel: +44 (0) 1403 799 190

Fax: +44 (0) 845 263 8907

medsales@gandn.com

www.gandn.com

